

Bilirubin Confirmation, Urine



Alternate Name:	Ictotest
Performing Lab:	New Hanover
Specimen Container:	Clean screw-cap container
Minimum Volume Required:	1.0 mL
Testing Availability	Routine: 24 hours/day Stat: No
Turnaround Time:	4 hours
Special Handling:	Transport specimen to the laboratory within 1 hour of collection. Avoid exposure to light.
Patient Preparation:	None
Specimen Stability:	
Reference Range:	Negative
Critical Value:	None
CPT Code:	81002
Testing Methodology:	Diazonium salt
Causes for Specimen Rejection:	Improper labeling Specimen too old Specimen exposed to light
Other Comments:	
Clinical Significance:	N/A