

Renal Panel

Acceptable Tube Types



Alternate Name:	None
Includes the following tests:	BUN Creatinine Sodium Potassium Chloride CO2 AGAP Calcium Phosphorus Glucose Albumin
Performing Lab:	New Hanover and Cape Fear
Specimen Container:	Yellow, green, or red top and microtainers
Minimum Volume Required:	1.0 mL
Testing Availability	Routine: 24 hours/day Stat: Yes
Turnaround Time:	Routine: 4 hours Stat: 1 hour
Special Handling:	None
Patient Preparation:	None
Specimen Stability:	N/A
Reference Range:	See individual tests
Critical Value:	See individual tests
CPT Code:	
Testing Methodology:	See individual tests
Causes for Specimen Rejection:	Improper labeling Hemolysis
Other Comments:	
Clinical Significance:	See individual tests

