

**Rubella IgG ( Immune Status )**



|                                       |                   |   |
|---------------------------------------|-------------------|---|
| <b>Alternate Name:</b>                |                   | Rubella Antibodies, German Measles Antibodies                                 |
| <b>Performing Lab:</b>                |                   | New Hanover   |
| <b>Specimen Container:</b>            |                   | Yellow serum separator tube or Red top tube or Tiger top Or Amber Microtainer |
| <b>Minimum Volume Required:</b>       |                   | 2 mL whole blood or 1 full microtainer  |
| <b>Testing Availability</b>           | <b>Routine:</b>   | Daily   |
|                                       | <b>Stat:</b>      | No  |
| <b>Turnaround Time:</b>               |                   | Within 24 hours   |
| <b>Special Handling:</b>              |                   | None  |
| <b>Patient Preparation:</b>           |                   | None  |
| <b>Specimen Stability:</b>            |                   |   |
| <b>Reference Range:</b>               |                   |   |
|                                       | <b>Negative:</b>  | < 5.0 IU/mL   |
|                                       | <b>Equivocal:</b> | > or equal to 5.0 IU/mL to 9.9 IU/mL  |
|                                       | <b>Positive:</b>  | > or equal to 10.0 IU/mL  |
| <b>Critical Value:</b>                |                   | None  |
| <b>CPT Code:</b>                      |                   | 86762   |
| <b>Testing Methodology:</b>           |                   | Chemiluminescence   |
| <b>Causes for Specimen Rejection:</b> |                   | Improper labeling   |
| <b>Other Comments:</b>                |                   |   |
| <b>Clinical Significance:</b>         |                   | Recommended for immune status determination                                   |